# dorsolite

REF.91111

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Back brace

Corset dorsal

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
Deutsch

**Please read the instructions carefully and thoroughly**

## **DECLARATION OF CONFORMITY**

As the manufacturer, ORTHOSERVICE AG declares with sole responsibility that this product is a Class I medical product and has been manufactured in accordance with the basic requirements of the directives 93/42/EEC and 2007/47/EEC. These operating instructions were created in accordance with the specification of the above-mentioned directives. They aim to ensure the safe and correct use of the medical product.

English

## **TRADEMARKS OF MATERIALS**

Velcro® is a registered trademark of Velcro® Industries B.V.;

## **SELECTION/SIZE**

| One size

| Length: 30 cm

## CHARACTERISTICS AND MATERIALS

- The support and straps are made from PE foam bonded with velveteen and a technical fabric.
- A central loop at the back provides a better grip when fastening.
- The straps can be adjusted and shortened using the Y-velcro straps:
  - The overall length can be adjusted;
  - The positioned of the shoulder straps can be adjusted.
- Fastens with a snap-fit buckle.
- Easy to put on, just like a backpack.

English

## INDICATIONS

- Helpful for correcting a hunchback, postural kyphosis and scapula winging.
- A useful aid for back pain caused by poor posture.

## CONTRAINDICATIONS

- No known


## WEARING THE BRACE FOR THE FIRST TIME

- 1 To put on the brace slide the patient's arms through the shoulder straps (fig A).
- 2 Adjust the position of the shoulder straps using the Y-tab on the top end (fig B):
  - Position the shoulder straps so that the padding protects the armpit.
  - If they are too long, remove the Y-Velcro® strap, cut it to the correct length and reposition the Y-Velcro® strap.
- 3 Fasten the snap-fit buckle at the front (fig C).
- 4 Adjust the length of the lower straps:
  - Stretch the straps to fit the patient (fig D).
  - If they are too long, remove the Y-Velcro® strap, cut it to the correct length and reposition the Y-Velcro® strap (fig H + fig I).
- 5 Place the brace back on the patient and adjust the fit if necessary.

## WEARING THE BRACE ONCE FITTED

- 6 Unfasten the snap-fit buckle.
- 7 To put on the brace slide the patient's arms through the shoulder straps (fig E).
- 8 Pull the straps at the sides (fig F) and fasten them at the front with the snap-fit buckle (fig G).

English

## MAINTENANCE

 do not bleach

 do not iron

 washing instructions: hand wash at max. 30° C with neutral soap.

Dry away from heat sources.

 no chemical cleaning

 do not tumble-dry

Do not dispose of the product or any of its components into the environment.

## PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. The device should not be worn in the vicinity of open flames.

## WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician and must be applied by an orthopaedic technician, the competent authority figure for the application as well as for information on the safe use, according to individual needs. To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. Use of the product is recommended for a single patient.

With hypersensitive patients, direct contact with skin may cause redness or irritation.

In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

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**Headquarter: ORTHOSERVICE AG**
**Niederlassung Deutschland: Orthoservice Deutschland GmbH**
Orthoservice AG (CH) · 6830 Chiasso (TI) · Switzerland