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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM**

Shoulder brace with immobilisation in external rotation at 15° or 30°

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician.

Please read the composition of the product on the inside label carefully. The device should not be worn in the vicinity of open flames or strong electromagnetic fields.

WARNINGS

The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician. The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements. The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. The manufacturer is no longer liable if the product is used or adapted inappropriately.

The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used.

SELECTION/SIZE

Item right	PR2-E113 15° · PR2-E115 30°		
Item left	PR2-E114 15° · PR2-E116 30°		
Size	S	M	L
Forearm length cm*	< 34	35/39	> 40
Colour	black		

*Measurement taken from elbow to end of closed fist
Specify right or left

MAINTENANCE

-  Do not bleach
-  No chemical cleaning
-  Do not iron
-  Do not tumble-dry
-  Washing instructions: Wash the cushion cover (after the foam inner padding has been removed), belts and the neck and underarm padding by hand with warm water (at max. 30°C) and neutral soap; rinse thoroughly. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Post-operative care (extra rotation) after operations on the rotator cuff, the joint of the glenoid labrum, capsular stabilisation, operations to the soft tissue of the shoulder
- Conservative treatment of initial stages of a luxation of the glenohumeral joint

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Cushion made from high-density foam. The fabric that comes into direct contact with the body has silver fibres woven into it
- Neck protection sleeve
- Covering and sleeve with a nylon gauze on the outside to facilitate the Velcro® grip
- Buckle straps and adjustable Velcro® fastening
- Option of extra rotation of 15° or 30°
- Rehabilitation ball
- Thumb support
- Arm support pocket made of breathable materi

PUTTING ON THE APPLIANCE

ADAPTING FOR THE DOCTOR/TECHNICIAN

- 1 Position the cushion on the side corresponding to the shoulder to be immobilised and place the arm in the appropriate breathing support pocket (fig. A). The arm-support pocket for the stabilisation of the limb may be repositioned on the cushion (using the Velcro®), if necessary.
- 2 Close the waistband with the large click-fasten buckle (fig. B) and tighten it with the corresponding Velcro® (fig. C). If it is too long, it may be shortened by cutting it to size and repositioning the Velcro® in a Y shape.
- 3 Close the arm-support pocket with the appropriate Velcro® strips (fig. D).
- 4 Position the thumb support (and the associated rehabilitation ball), using the Velcro® to attach it to the cushion (fig. E).
- 5 For correct extra-rotation positioning, position yourself in front of a table and make sure that the white line on the upper part of the cushion is parallel to the table.
- 6 Slip the opposite arm into the triangle created by the armpit and the clavicle strap, so that the orange-marked strap is above the arm and the white-marked one passes below the armpit (fig. F).
- 7 Close the guide by attaching the small click-fasten buckle (fig. G).
- 8 Tighten the (orange) shoulder strap with the appropriate Velcro® (fig. H). If it is too long, it may be shortened by cutting it to size and repositioning the Velcro® in a Y shape.
- 9 Tighten the (white) under-armpit strap with the appropriate Velcro® (fig. I). If it is too long, it may be shortened by cutting it to size and repositioning the Velcro® in a Y shape.
- 10 Check and if necessary reposition the sleeve padding at the level of the neck and under the armpit.
- 11 Tighten the rear strap with the appropriate Velcro® (fig. K). If it is too long, it may be shortened by cutting it to size and repositioning the Velcro® in a Y shape.

PUTTING ON THE APPLIANCE FOR THE PATIENT

Maintain the configuration set when using for the first time. The appliance may be removed by opening the arm-support pocket and unfastening the closure buckles. The support is reapplied as follows:

- 1 Position the cushion on the side corresponding to the shoulder to be immobilised and place the arm in the appropriate transpiring support pocket (fig. A).
- 2 Close the waistband with the large click-fasten buckle (fig. B).
- 3 Close the arm-support pocket with the appropriate Velcro® strips (fig. D) and slip the thumb into the appropriate support.
- 4 For correct extra-rotation positioning, position yourself in front of a table and make sure that the white line on the upper part of the cushion is parallel to the table.
- 5 Slip the opposite arm into the triangle created by the armpit and the clavicle strap, so that the orange marked strap is above the arm and the white-marked one passes below the armpit (fig. F).
- 6 Close the guide by attaching the small click-fasten buckle (fig. G).

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Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.