
CENT

ORTHOSERVICE

RO+TEN Take care feel better

OBJEC

REF. 460

REF. 463

REF.404N

REF. 416

Short tubular knee support

Neoprene short tubular knee

Neoprene tubular knee support

support with open kneecap

with patella stabiliser

Knee patella stabiliser

Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89 info@orthoservice.com · www.orthoservice.com Niederlassung Deutschland: Orthoservice Deutschland GmbH

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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM

REF. 460 Short tubular knee support

REF. 463 Neoprene short tubular knee support with open kneecap

REF. 404N Neoprene tubular knee support with patella stabiliser

REF. 416 Knee patella stabiliser

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class medical device that has been manufactured in accordance with the requirements of Directives 93/42/ EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

Velcro® is a registered trademark of Velcro Industries B.V.

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. If in doubt about how to apply, contact an orthopaedic technician. The device should not be worn in the vicinity of open flames. People prone to dermatitis or allergies to products resulting from the neoprene must not use this brace.

The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician. The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician.

Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. The manufacturer is no longer liable if the product is used or adapted inappropriately. The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals.

In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used. SELECTION/SIZE

Code	REF. 460 - REF. 463							
Size	XS	S	M	L	XL	XXL		
Knee center circumf. cm	27/30	31/33	34/36	37/39	40/43	44/47		
Lenght cm	26	26	26	26	26	26		
Colour	black REF. 460 - black and blue REF. 463							

Fits right and left

Code	REF. 404N - REF. 416							
Size	XS	S	M	L	XL	XXL		
Knee center circumf. cm	30/33	33/36	36/39	39/43	43/48	48/53		
Lenght cm	33	33	33	33	33	33		
Colour	black and blue							
						i		

Fits right and left

MAINTENANCE

⋈ No chemical cleaning

■ Do not tumble-dry

Washing instructions: Remove the palm splint from its seat. Hand wash in warm water with mild soap: rinse it thoroughly. Allow to dry away from heat sources

Do not dispose of the product or any of its components into the environmen

INDICATIONS REF. 460/463/404N

- Knee pain caused by bruising
- Early stage arthrosis
 REF. 460/463 only Grade 1 collateral ligament sprains
- REF. 460 only: Tendonitis
 REF. 463 only: Prevention of patellar dislocation during sports activities
- REF. 404N only
- · Femoro-patellar chondropathy and chondromalacia in the knee
- · Mild instability of the patella
- · Prevention of displacement of the patella during sporting activity

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Currently no known

FEATURES AND MATERIALS

- Double-lined orthonaedic neonrene
- Lined with soft velour for enhanced comfort and to help avoid excessive sweating
- · REF. 463 only: Patellar opening with reinforced edging for patellar stabilisation
- REF. 404N only
- · A pocket around the knee for positioning the horseshoe stabiliser
- REF. 416 only · Open kneecap with a lateral J-shaped foam stabiliser
- Straps with Velcro® fasteners
- · Medial and lateral spiral stavs

PUTTING ON THE APPLIANCE REF. 460/463

- 1 Slide the tubular support onto the knee with the writing toward the top.
 2 **Model REF. 460** (fig. A): Position the knee support so that the seams are at the side.
- Model REF. 463 (fig. B): Position the knee support so that the kneecap is perfectly centred in the

opening provided.
PUTTING ON THE APPLIANCE REF. 404N

1 Slide the tubular support onto the knee with the writing toward the top and position the knee support so that the kneecap is perfectly centred in the opening provided.

2 If necessary, position the horseshoe stabiliser around the kneecap in the pocket provided

2 Slide the tubular support onto the knee (fig. D) and position it so that the kneecap is perfectly centred

PUTTING ON THE APPLIANCE REF. 416 1 Undo the patella straps.

in the opening provided and the stabiliser is to the side of the kneecap (fig. E).

3 Pull the patella straps and fasten them with Velcro® (fig. FGD). N.B.: if the Velcro® is too long and pokes out, it can be cut to size

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