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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM

Night positioning ankle brace

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. Please read the composition of the product on the inside label carefully. The device should not be worn in the vicinity of open flames or strong electromagnetic fields.

WARNINGS

The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician. The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements. The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function.

Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. The manufacturer is no longer liable if the product is used or adapted inappropriately.

The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used.

SELECTION/SIZE

Code	REF.90062	
Size	S/M	L/XL
Shoe size	36/42	43/48
Right and left		

MAINTENANCE

- Do not bleach
- No chemical cleaning
- Do not iron
- Do not tumble-dry
- Washing instructions: hand wash at max. 30° C with neutral soap.
Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Plantar fascitis
- Achilles tendinopathy
- Foot drop
- Pain upon awakening

CONTRAINdications

Currently no known

CHARACTERISTICS AND MATERIALS

- The hinged polymer splint can be positioned for personalised support using the adjustable plantar flexion strap
- Comfortable, soft padding
- The adjustable Velcro® calf and foot straps are padded for optimum comfort.

PUTTING ON THE APPLIANCE

- 1 Undo all the straps and fasten them back on themselves for the moment.
- 2 Slide the foot into position and rest the splint of the top of the foot and the ankle (fig. A).
- 3 Pull the Velcro® calf strap and fasten (fig. B).
- 4 Next pull the foot (fig. C) strap and fasten.
- 5 [For the first time of wearing only] Adjust the dorsal flexion using the Velcro® strap (fig. D).
- 6 If greater stabilisation is required, an additional strap is provided that can be placed around the ankle (fig. E).

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