RADIO LUCENT

LOW PROFILE

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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM

Ankle foot orthosis

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/EEC and 2007/47/CE. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. Please read the composition of the product on the inside label carefully. The device should not be worn in the vicinity of open flames or strong electromagnetic fields.

WARNINGS

The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician. The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements. The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be the last word. The manufacturer is no longer liable if the product is used or adapted inappropriately. The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used.

SELECTION/SIZE

Code	REF.9050 - REF.9051
Size	S M L XL
Height cm	28 29 34 39
Shoe size	35/39 40/43 41/44 >44
Colour	REF.9050 = white - REF.9051 = light blue

specify right or left

MAINTENANCE

- Do not bleach
- Do not iron
- Do not tumble-dry
- Washing instructions: hand wash at max. 30 °C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources (or, better, dry with a cloth)
- Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Weakness of dorsiflexion of the foot (drop foot)
- Light to moderate spastic equinus foot

CONTRAINdications

Currently known

CHARACTERISTICS AND MATERIALS

- Body REF.9050 of polypropylene structure with rounded-off edges
- Body REF.9051 made of hydroscopic, non-malleable, robust and elastic polyoxymethylene (acetal resin)
- Soft padding for protection of the tibia and the calf (MTP padding, with light velour covering on both sides)
- Velcro® closure at the calf
- Thermo-formable (REF.9050)
- Velcro® closure across the arch of the foot

ADAPTING / PUTTING ON THE APPLIANCE

ADAPTING FOR THE DOCTOR / TECHNICIAN

- Where necessary, trim the end of the foot section (fig. A) with scissors, so as to adapt the appliance to the shoe and foot of the patient. Roll over the edges.
- REF.9050: If necessary, thermoform the appliance, in order to ensure that it fits the patient's leg better.
- Insert the appliance into the shoe.
- Put the shoe on with the appliance already inserted.
- Fasten the padding on the calf (fig. B).
- Pull the Velcro® strap around the calf (fig. C).
- Where necessary, use the foot strap to secure the appliance (fig. D).

PUTTING ON THE APPLIANCE FOR THE PATIENT

- Insert the appliance into the shoe.
- Put the shoe on with the appliance already inserted.
- Fasten the padding on the calf (fig. B).
- Pull the Velcro® strap around the calf (fig. C).
- Where used, fasten the foot strap (fig. D).

PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM

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